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Usage Instruction

為了您的安全，請您務必花費幾分鐘的時間，在充分了解本說明書之後，按照內容作正確的使用。
In order to protect your safety, please spend some times to read thought this instruction and operate the device with the direction before use.

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4 712702 487031 0

KU.KU Duckbill Storage Milk Bottle Set- 3 Bottles

When mothers have collected excess breast milk, we suggest to change a new storage milk bottle to collect breast milk in order to manage the amount of breast milk easily. KU.KU Duckbill storage milk bottle set- 3 bottles are designed for mothers to storage more breast milk.

- The milk bottle can use standard nipple.
- The milk bottle can be placed on the breast pump to used as the storage milk bottle.
- Put on the nipple and cover to be used as the feeding bottle.
- Foodstuff level PP material can be sterilized by high temperature.
- It can be repeat for use which is money saving and environment friendly.

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No.27, Lane 290, Sec.2, Jia Dong Road, Chang Hua City 500-84, Taiwan

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Suggestions:

At room temperature-for 8 hours at a room temperature of 25 °C
In the refrigerator for up to 5 days at a refrigerator temperature of 0 °C to 4 °C
In the freezer for up to 2 weeks in a freezer compartment contained inside a refrigerator, for 3 to 6months inside a freezer compartment in a refrigerator with a separate freezer door, and 6 months or longer if stored in a separate deep freezer. Do not store breast milk in the freezer door rack. Milk stored in the door is exposed to a wide temperature changes due to the frequent opening and closing of the freezer door.

◎ 儲乳使用說明 / Instruction of the storage milk bottle

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1. After finishing the breast milk collection, release the milk bottle.
2. Place the milk bottle then it can be use as the storage milk bottle.
3. If the storage milk bottle is used for feeding, please remove the milk storage spacer then put on the pacifier and milk bottle cover then it can be used for feeding.
4. The storage milk bottle is used to storage the excess breast milk in the fridge (use within 24 hours) and warns up by hot water to feed baby when it is needed. (Test the temperature before feeding.)
5. Ku Ku duckbill milk bottle are suggested to be used. It can be fit to the pump better, avoid the risk of fall and easy to storage.

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Breast-Feeding Best for Babies

The breast milk is the best food for baby even the humanized milk is not able to compare with. Baby feed by breast milk is usually grow more health then the bottle feeding baby therefore smart Mon shall chose breast-feeding.. Now born baby has three basic needs: warn embrace from Mon, the breast milk from Mon and the sense of security from Mon's company that can be done by every single mother.

The benefits from breast feeding to baby

1. Breast milk is the most nature and accord with human physical baby food.
2. Breast milk is fresh and in stable temperature.
3. Breast milk is able to reduce the symptoms of baby's allergy such as ecaema, nettle rush, asthma, allergic rhinitis and etc.
4. The possibility is lower for breast feeding baby to get cold, diarrhea, constipation, skin problem and etc then the bottle feeding one. This is due to the breast milk contained various antibody to enhance the resistance which is especially important to the baby.
5. Breast milk is clear and easy to digest which can reduce the distention and uncomfortable.
6. There are less malnutrition and over nutrition on breast feeding baby.

The important colostrum

Colostrum is the earliest breast milk lactate by mother. The amount of the colostrum is very little. The light yellow color milk contains very rich of nutrition and antibody which can help baby on digestion.

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1. Power save

KU.KU Duckbill manual breast pump with handle which save power to operate it. Easy to use and assemble. It is ideal for reducing uncomfortable for the breast feeding mother.

2. Smart suction and release

The suction can according to the mother needs to control the movement, let the mother feel very comfortable.

3. Comfortable massage

The massage cushion with embossed petal air room design that massage the area behind of nipple to stimulate the milk flow.

4. Easy to use


The mother only uses one hand to operate and avoid hand fatigue, whole design meets Human Ergonomics, curved handle for either hand to operate it.

5. Low sound and quiet

The baby won't interrupt during processing of breast pump because it is low sound.

6. Health protection


The design of breast shield cap is convenient to carry and keep the shield clean.

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
- Clean both hands and breast before use.
- Breast pump must be well sterilized and assembled.
- Fit the milk bottle on the breast pump for collecting breast milk.
- The manual handle can be adjusted to the best angle for mother by turning by left hand (picture 1) and right hand (picture 2).
- Place the breast pump on the breast properly and press the handle on the pump.
- Press the handle briskly to massage the breast, press to half and release until the breast milk come out.
- Slow down the speed after the breast milk start to come out. Press down the handle and imitate the nature temple of baby's suction to keep the press and release in the same frequency.
- When it is not able to collect the milk, please Check is the pump cover, plate, suction ring and suction ring base fit properly.

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◎ 手動組裝說明 / Assembling instruction of manual breast pump

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- Place the breast pump plate ② into breast pump cover ③ and make it completely close together.
- Place the suction ring ⑥ into the silica base ⑤ from bottom to the top.
- Place the base ⑦ on the breast pump cover and turn the breast pump cover to go with the tenon A.
- After the set up, place the suction ring base ⑤ ⑥ on the manual appliance base ⑦.
- Place the manual handle front part ⑪ on the suction ring and fix the handle on the appliance base and put on the handle cover ④.
- Place the anti-back flow silica tube ⑨ and milk bottle gasket ⑧ and close the milk bottle and breast pump cover.

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- Before first use and every use must wash all the parts with warm water. The pacifier, milk bottle, breast pump cover, cover plate, anti-back flow silicate tube, milk bottle spacer and so on are sterilized.
- Please clean and sterilize the pump before using Clean it up immediately after using and keep the container dry.
- Please do not expose this product under sunshine.
- Please follow up the assembly or disassembly order completely. Please do not lose small parts when you clean it.